

# skin peeler by REVIDERM – the original unique design meets modern technology

Microdermabrasion is the core of the regeneration treatments from REVIDERM skintelligence. Originally developed for the medical field, today it is the leading non-operative anti-aging method found in the professional beauty sector.

The success of controlled microdermabrasion lies in the functionality of the REVIDERM Skin Peeler. The upper skin layers are gently removed by a micro-fine stream of crystals in a controlled way. Your customers will be convinced even after the first treatment by their soft and balanced skin.

At the same time the skin is ideally prepared for the penetration of active agents. Valuable essences from the REVIDERM skintelligence system guarantee an optimal active agent supply into the skin. The combination of optimal peeling and vacuum activates further biological repair processes in the skin. The skin's own collagen production is stimulated and microcirculation is increased. This leads to an increased production of young, functional skin cells in the basal membrane. These migrate to the upper layers and continuously care for the regeneration of the skin from within.


As leading worldwide manufacturer REVIDERM sets technological standards. Numerous international commendations confirm the highest quality. Our know-how and our service-oriented approach guarantee you absolute customer satisfaction and reliability.

Microdermabrasion is outstandingly suitable for all regeneration treatments. The broad spectrum of indications makes these treatment methods indispensable for your beauty clinic and for enthralled customers.


- horny cells (hyperkeratosis)
- atrophic skin
- light and sun-damaged skin
- pigmentation disorders and age spots

- impurities and acne scars
- body treatments striae, stretch marks and cellulite
- · scars on the face and body
- pre and post treatment for facial injections


For perfect results microdermabrasion can be excellently combined with ultrasound. This increases the effects of microdermabrasion and the skin is rejuvenated in a natural way.



 Enzymes split the dead skin cells and prepare the skin for microdermabrasion



 Microdermabrasion abrodes the dead skin cells gently and safely. The skin surface appears smooth and tender



 Individually chosen active agents are brought deeply into the skin, where they unfold their maximum effectiveness.

## REVIDERM innovation guarantee

All REVIDERM Skin Peelers of the new generation are based on the latest technology and guarantee

- · gentle and effective peeling
- safe and controlled abrasion
- · quiet and relaxed treatment
- dust-free and clean handling





### The Entry Model

practical characteristics for professional handling

#### individual functions

- electronic vacuum regulation (4 adjustable levels)
- · diagnostic mode of the vacuum status
- · tube-in-tube system
- no dust disposable crystal containers



#### The Professional Model

consequent advancement of professional handling

#### individual intensity

- electronic vacuum regulation (8 adjustable levels)
- crystal regulation (4 adjustable levels)
- · diagnostic mode to test the vacuum and filter status
- · automatic filter cleaning
- · humidity filter system to prevent the penetration of moisture
- tube-in-tube system
- no dust disposable crystal containers



#### The Ultimate Model

all-in-one device combining microdermabrasion and ultrasound

#### as the Skin Peeler New Generation including ultrasound unit

- with 2 ultrasound frequencies: 1 MHz for facial treatments and 3 MHz for body treatments
- with pulsed and permanent sound adjustable in 3 different levels
- · timer function to pre-set the individual treatment duration
- · titanium ultrasound hand-pieces for body and face

Exclusively Distributed by



Dubai

+971 4 3884549

Abu Dhabi

+971 2 6275926

Oatar

+974 33921190

www.specbritegroup.com

ISO 9001:2015 CERTIFIED (CRS) ((1)







